

CHLORAXiD 2%

Liquid for root canals rinsing

CE₂₄₀₉

PRODUCT FOR DENTAL USE ONLY

Package:

Bottle of 200g, 400g

Composition: sodium hypochlorite, water Ph. Eur.

Intended use

During mechanical canal widening it removes the non-vital pulp debris.

It cleans the canal, removes the smear layer in order to expose the dentin root canal orifices.

It prevents from the tooth discoloration that might occur after filling the non-rinsed root canal.

Adapter with a plug for taking the liquid preparations into the syringe

1. Unscrew the bottle cap

2. Place adapter into the bottle neck and press it firmly.

3. In order to use the bottle remove the blue plug, screw the end of the syringe into the adapter's opening, slope the bottle with preparation upside down and take the liquid into the syringe, while holding down the adapter at the same time. Next put the bottle upright on the table top and gently unscrew the filled syringe from adapter. It is recommended to use luer-lock syringes with rubber plunger.

4. After each use close the adapter's opening with a blue plug.

5. For transportation remove adapter and twist the bottle with a bottle cap.

Instruction for use

Insert the liquid CHLORAXID 2% into the root canal by means of the syringe with needle with side portal.

Rinse the canal copiously, approximately 10 ml for one canal.

Sodium hypochlorite contained in the product inactivates immediately in contact with organic substance.

While concentration of sodium hypochlorite is decreasing, time of rinsing and volume should be necessarily simultaneously increased.

Cautions

During root canal preparation do not use too heavy stream of the liquid because it can cause the liquid penetration to the periapical tissues and damage them, characterized by pain and swelling of soft tissues.

NOTICE!!! CHLORAXID 2% can irritate soft tissues of the oral cavity, eyes and skin. Necessary use of rubber-dam while working with the product.

Recommended use of personal protective equipment (glasses, protective visor, gloves). Avoid contact with eyes, soft tissues of oral cavity and skin.

In case of accidental contact rinse with plenty of water.

Contra-indications

Do not use product CHLORAXID 2% for patients with hypersensitivity to the product's components.

Storage

Store in the original containers at temperature 2 – 8 °C (in the refrigerator). Store in a securely closed packaging. Protect the product from the sunlight. Keep out of reach of children.


Transportation

The product manufacturer is admitting the product transportation within the period not longer than 10 days at temperature not exceeding 40°C. Short transitional period at given temperature range has no influence on safety and effectiveness of the product throughout its shelf life.

Do not use the product after the expiration date placed on the container.

Once container is open, the expiry date doesn't change if it is securely closed after each use.

Empty containers should be disposed of responsibly.



PPH CERKAMED Wojciech Pawłowski
37-450 Stalowa Wola, Polska
ul. Kwiatkowskiego 1
Phone/fax 0048 15 842-35-85