PL

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
EN

Biomaterial for filling and rebuilding of root canals

**PLEASE READ THE FOLLOWING INSTRUCTIONS BEFORE USE PRODUCT FOR DENTAL USE ONLY**

#### INGREDIENTS:

*Bio MTA+ powder*: calcium oxide, hydroxyapatite, oxides of: silicon, iron, aluminium, sodium, potassium, bismuth, magnesium, zirconium; calcium phosphate.  
*Bio MTA+ liquid*: Ph Eur purified water, calcium catalyst.

#### INDICATIONS FOR USE

The product is intended for dental treatment as a material for filling and remineralization of root canals. When mixed with the *Bio MTA+ liquid*, the working life is ca. 4 min, and the compound sets completely in ca. 2 h. The high content of calcium ions remineralises the tooth tissue. The silicon and calcium compounds support tissue regeneration after wall perforation and intracanal resorption. *Bio MTA+* contains hydroxyapatite, a natural component of bones, which integrates the compound perfectly into the bone structure.  
*Bio MTA+* is especially recommended in:  
- root canal wall perforation  
- intracanal resorption  
- filling of the posterior root apex  
- direct coating of the pulp  
- pulp amputation  
- treatment of teeth with interrupted root development  
*Bio MTA+* is not suitable as a standalone and permanent root canal filling. Once its two components have been mixed, *Bio MTA+* exhibits extremely plastic and cohesive consistence without brittling and stays on the instruments. This facilitates the application *Bio MTA+* when compared with the traditional MTA compounds.

#### INSTRUCTIONS FOR USE

Cut off the *Bio MTA+ liquid* pipette tip. Move the entire contents of the *Bio MTA+ powder* glass vial with 1 to 2 drops of the *Bio MTA+ liquid* on a mixing plate. Mix for ca. 30 s until the compound reaches a consistence of soft plasticine (modelling putty). If the compound mix is too thick or too brittle, add one more drop of the *Bio MTA+ liquid*.  
**CAUTION!** Mix the *Bio MTA+ powder* with the *Bio MTA+ liquid* only. Do not replace the liquid with water or other fluids. Place the compound at the target site with an applicator. One container of the *Bio MTA+ powder* contains a one-off dose. Once mixed with the *Bio MTA+ liquid*, consume the mixed dose within a maximum of 4 minutes. The compound cures and hardens afterwards.

#### CONTRAINDICATIONS

Do not use in patients who are hypersensitive to any of the ingredients. Do not use in inflammatory conditions. Acidic pH values disrupt the setting of *Bio MTA+*. Do not mix the *Bio MTA+ powder* with water or other fluids. Use the *Bio MTA+ liquid* only.

#### PRECAUTIONS

The product can cause eye damage. The product irritates the oral mucosa, the skin and eyes. In case of contact, rinse with plenty of water and contact a doctor immediately. If swallowed, do not induce vomiting. Drink plenty of water. Contact a doctor immediately. Use rubber dams, disposable gloves and safety goggles when working with the product.

#### STORAGE


Store in a tightly closed original container below 25 °C. Protect from moisture. Keep out of the reach of children. The expiry date is on the direct packaging.

#### DISPOSAL OF PACKAGING

Waste containers should be disposed of or returned to the manufacturer.

#### PACKAGING

Glass vial containing 0.14 g of the *Bio MTA+ powder*; pipette containing 1 ml of the *Bio MTA+ liquid*.

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